



# **TissueMAPS Toolbox Documentation**

*Release 0.1.0*

**Markus D. Herrmann, Robin Hafen**

June 27, 2015



## CONTENTS

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Subpackages . . . . .	3
1.2	Configuration settings . . . . .	3
1.3	Documentation . . . . .	4
<b>2</b>	<b>Indices and tables</b>	<b>5</b>



Contents:



## INTRODUCTION

**tmt** is a Python package that bundles image processing and data analysis tools for [TissueMAPS](#).

### 1.1 Subpackages

**align** - Alignment of images between different acquisition cycles.

**corilla** - Correction of illumination artifacts.

**dafu** - Data fusion for [Jterator](#) projects.

**illuminati** - Creation of image pyramids.

**visi** - Conversion of Visitron's STK files to PNG images with optional renaming.

### 1.2 Configuration settings

Configurations are defined in `.config` [YAML](#) files to specify the experiment layout, such as the directory structure on disk.

Paths and filenames are described with [Python format strings](#). The **replacement fields** surrounded by curly braces `{ }` are then automatically replaced with experiment specific variables.

**To this end, you can use the following *replacement fields*:**

- *experiment\_dir*: absolute path to the experiment directory
- *experiment*: name of the experiment folder
- *subexperiment*: name of a subexperiment folder, i.e. a subfolder of the experiment folder
- *cycle*: number of a subexperiment
- *channel*: number of a channel of intensity images (layers)
- *objects*: name of objects in segmentation images (masks)

```
SUBEXPERIMENTS_EXIST: Yes

# Path format strings
IMAGE_FOLDER_LOCATION: '{experiment_dir}/{subexperiment}/images'
SHIFT_FOLDER_LOCATION: '{experiment_dir}/{subexperiment}/shift'
STATS_FOLDER_LOCATION: '{experiment_dir}/{subexperiment}/stats'
SEGMENTATION_FOLDER_LOCATION: '{experiment_dir}/{subexperiment}/segmentations'
LAYERS_FOLDER_LOCATION: '{experiment_dir}/layers'
DATA_FILE_LOCATION: '{experiment_dir}/data.h5'
```

```
# Filename format strings
SUBEXPERIMENT_FOLDER_FORMAT: '{experiment}_{cycle:0>2}'
SUBEXPERIMENT_FILE_FORMAT: '{experiment}_{cycle}'
STATS_FILE_FORMAT: 'illumstats_channel{channel:0>3}.h5'
SHIFT_FILE_FORMAT: 'shiftDescriptor.json'

# Regular expression patterns to extract information encoded in filenames
EXPERIMENT_FROM_FILENAME: '^([\_]+)'
CYCLE_FROM_FILENAME: '_(\d+)_ '
COORDINATES_FROM_FILENAME: '_r(\d+)_c(\d+)_ '
COORDINATES_IN_FILENAME_ONE_BASED: Yes
SITE_FROM_FILENAME: '_s(\d+)_ '
CHANNEL_FROM_FILENAME: 'C(\d+)\.png$'
OBJECTS_FROM_FILENAME: '_segmented(\w+)\.png$'

# Should Vips image processing library be used? Required for pyramid creation!
USE_VIPS_LIBRARY: Yes

# These settings are hard-coded in TissueMAPS, so don't change them!
LAYERS_FOLDER_LOCATION: '{experiment_dir}/layers'
ID_TABLES_FOLDER_LOCATION: '{experiment_dir}/id_tables'
ID_PYRAMIDS_FOLDER_LOCATION: '{experiment_dir}/id_pyramids'
DATA_FILE_LOCATION: '{experiment_dir}/data.h5'
```

NOTE: Quotes are generally not required around strings in YAML syntax, but are necessary here because of the curly braces in the format strings!

## 1.3 Documentation

Sphinx is used for the documentation of source code in combination with the [Napoleon extension](#) to support the [reStructuredText NumPy style](#).

Documentation is located under *docs* and will ultimately be hosted on [Read the Docs](#).

To update the documentation upon changes in the source code, do

```
sphinx-apidoc -o ./docs ./tmt
```

To build HTML, do

```
cd docs
make html
```

The generated HTML files are located at *docs/\_build/html*.



## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`